

# CSH1 Antibody



PACO08627

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CSH1

**Uniprot**

P01243

**Synonyms:**

chorionic somatomammotropin hormone 1 (placental lactogen); CSH1; CSA; CSMT; FLJ75407; PL ;

**Immunogen:**

Human CSH1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A