

CSPG5 Antibody



PACO08636

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CSPG5

Uniprot

O95196

Synonyms:

chondroitin sulfate proteoglycan 5 (neuroglycan C); CSPG5; MGC44034; NGC ;

Immunogen:

Human CSPG5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A