

CST8 Antibody



PACO08644

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CST8

Uniprot

O60676

Synonyms:

cystatin 8 (cystatin-related epididymal specific); CST8; CRES ;

Immunogen:

Human CST8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A