

CTAG1A Antibody



PACO08648

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CTAG1A

Uniprot

P78358

Synonyms:

cancer/testis antigen 1A; CTAG1A; ESO1; LAGE2A ;

Immunogen:

Human CTAG1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A