

CTAGE1 Antibody



PACO08649

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CTAGE1

Uniprot

Q96RT6

Synonyms:

cutaneous T-cell lymphoma-associated antigen 1; CTAGE1; CTAGE-1; CTAGE-2 ; CTAGE; CT21.1; CT21.2, CTAGE2

Immunogen:

Human CTAGE1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A