

PACO08665

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CTPS2

Uniprot

Q9NRF8

Synonyms:

CTP synthase II; CTPS2; DKFZp686C17207; FLJ43358; MGC32997 ;

Immunogen:

Human CTPS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A