

CTSH Antibody



PACO08672

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CTSH

Uniprot

P09668

Synonyms:

cathepsin H; CTSH; ACC-4; ACC-5; CPSB; DKFZp686B24257; MGC1519; minichain ;

Immunogen:

Human CTSH.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A