

PACO08673

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CTSK

**Uniprot**

P43235

**Synonyms:**

cathepsin K; CTSK; CTS02; CTSO; CTSO1; CTSO2; MGC23107; PKND; PYCD ;

**Immunogen:**

Human CTSK.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A