

CTTNBP2 Antibody



PACO08676

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CTTNBP2

Uniprot

Q8WZ74

Synonyms:

cortactin binding protein 2; CTTNBP2; C7orf8; CORTBP2; FLJ34229; KIAA1758; MGC104579; Orf4 ;

Immunogen:

Human CTTNBP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A