

# CXCL1 Antibody



PACO08685

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CXCL1

**Uniprot**

P09341

**Synonyms:**

chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha); CXCL1; FSP; GRO1; GROa; MGSA; MGSA-a; NAP-3; SCYB1 ;

**Immunogen:**

Human CXCL1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A