

CXCL14 Antibody



PACO08688

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CXCL14

Uniprot

O95715

Synonyms:

chemokine (C-X-C motif) ligand 14; CXCL14; BMAC; BRAK; KS1; Kec; MGC10687; MIP-2g; NJAC; SCYB14; bolekin ;

Immunogen:

Human CXCL14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A