

# CXCL5 Antibody



PACO08691

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CXCL5

**Uniprot**

P42830

**Synonyms:**

12kd(8 or 8.3kd);

**Immunogen:**

Human CXCL5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A