

CYFIP2 Antibody



PACO08702

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CYFIP2

Uniprot

Q96F07

Synonyms:

cytoplasmic FMR1 interacting protein 2; CYFIP2; PIR121 ;

Immunogen:

Human CYFIP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A