

CYLC2 Antibody



PACO08704

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CYLC2

Uniprot

Q14093

Synonyms:

cylicin, basic protein of sperm head cytoskeleton 2; CYLC2; MGC129591 ;

Immunogen:

Human CYLC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A