

# CYP26B1 Antibody



PACO08711

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

CYP26B1

**Uniprot**

Q9NR63

**Synonyms:**

cytochrome P450, family 26, subfamily B, polypeptide 1; CYP26B1; CYP26A2; DKFZp686G0638; MGC129613; P450RAI-2 ;

**Immunogen:**

Human CYP26B1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A