

# CYP2A6 Antibody



PACO08713

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CYP2A6

**Uniprot**

P11509

**Synonyms:**

cytochrome P450, family 2, subfamily A, polypeptide 6; CYP2A6; CPA6; CYP2A; CYP2A3; P450C2A; P450PB ;

**Immunogen:**

Human CYP2A6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A