

# CYP2F1 Antibody



PACO08715

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CYP2F1

**Uniprot**

P24903

**Synonyms:**

cytochrome P450, family 2, subfamily F, polypeptide 1; CYP2F1; C2F1; CYP2F;  
MGC126121 ;

**Immunogen:**

Human CYP2F1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A