

CYP2W1 Antibody



PACO08718

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CYP2W1

Uniprot

Q8TAV3

Synonyms:

cytochrome P450, family 2, subfamily W, polypeptide 1; CYP2W1; MGC34287 ;

Immunogen:

Human CYP2W1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A