

CYP3A5 Antibody



PACO08720

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CYP3A5

Uniprot

P20815

Synonyms:

cytochrome P450, family 3, subfamily A, polypeptide 5; CYP3A5; CP35; P450PCN3; PCN3 ;

Immunogen:

Human CYP3A5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A