

CYTH2 Antibody



PACO08728

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CYTH2

Uniprot

Q99418

Synonyms:

cytohesin 2; CYTH2; ARNO; CTS18; CTS18.1; PSCD2; PSCD2L; SEC7L; Sec7p-L; Sec7p-like ;

Immunogen:

Human CYTH2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A