

DAD1 Antibody



PACO08735

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

DAD1

Uniprot

P61803

Synonyms:

12kd;

Immunogen:

Human DAD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A