

DAG1 Antibody



PACO08736

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DAG1

Uniprot

Q14118

Synonyms:

dystroglycan 1 (dystrophin-associated glycoprotein 1); DAG1;156DAG; A3a; AGRNR; DAG ;

Immunogen:

Human DAG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A