

# DCLRE1C Antibody



PACO08759

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

DCLRE1C

**Uniprot**

Q96SD1

**Synonyms:**

DNA cross-link repair 1C (PSO2 homolog, *S. cerevisiae*); DCLRE1C; A-SCID; DCLREC1C; FLJ11360; FLJ36438; RS-SCID; SCIDA; SNM1C ;

**Immunogen:**

Human DCLRE1C.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A