

# DCN Antibody



PACO08760

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DCN

**Uniprot**

P07585

**Synonyms:**

decorin; DCN; CSCD; DSPG2; PG40; PGII; PGS2; SLRR1B ;

**Immunogen:**

Human DCN.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A