

DCUN1D5 Antibody



PACO08769

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

DCUN1D5

Uniprot

Q9BTE7

Synonyms:

DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae); DCUN1D5; FLJ32431; MGC2714 ;

Immunogen:

Human DCUN1D5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A