

DDB1 Antibody



PACO08774

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DDB1

Uniprot

Q16531

Synonyms:

damage-specific DNA binding protein 1, 127kDa; DDB1; DDBA; UV-DDB1; XAP1; XPCE; XPE; XPE-BF ;

Immunogen:

Human DDB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A