

DDIT4L Antibody



PACO08779

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DDIT4L

Uniprot

Q96D03

Synonyms:

DNA-damage-inducible transcript 4-like; DDIT4L; REDD2; Rtp801L ;

Immunogen:

Human DDIT4L.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A