

DDX19B Antibody



PACO08783

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DDX19B

Uniprot

Q9UMR2

Synonyms:

DEAD (Asp-Glu-Ala-As) box polypeptide 19B; DDX19B; DBP5; DDX19; RNAh ;

Immunogen:

Human DDX19B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A