

DDX27 Antibody



PACO08789

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DDX27

Uniprot

Q96GQ7

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 27; DDX27; DKFZp667N057; FLJ12917; FLJ20596; FLJ22238; HSPC259; MGC1018; MGC163147; PP3241; RHL P; Rrp3p; dJ686N3.1 ;

Immunogen:

Human DDX27.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A