

DDX43 Antibody



PACO08793

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DDX43

Uniprot

Q9NXZ2

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 43; DDX43; CT13; DKFZp434H2114; HAGE ;

Immunogen:

Human DDX43.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A