

DDX50 Antibody



PACO08795

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DDX50

Uniprot

Q9BQ39

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 50; DDX50; GU2; GUB; MGC3199; RH-II/GuB ;

Immunogen:

Human DDX50.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A