

DEAF1 Antibody



PACO08799

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DEAF1

Uniprot

O75398

Synonyms:

deformed epidermal autoregulatory factor 1 (Drosophila); DEAF1; NUDR; SPN; ZMYND5 ;

Immunogen:

Human DEAF1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A