

DES Antibody



PACO08807

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DES

Uniprot

P17661

Synonyms:

desmin; DES; CMD1I; CSM1; CSM2; FLJ12025; FLJ39719; FLJ41013; FLJ41793 ;

Immunogen:

Human DES.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A