

# DGCR2 Antibody



PACO08810

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DGCR2

**Uniprot**

P98153

**Synonyms:**

DiGeorge syndrome critical region gene 2; DGCR2; DGS-C; DKFZp686I1730; IDD; KIAA0163; LAN; SEZ-12 ;

**Immunogen:**

Human DGCR2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A