

DGUOK Antibody



PACO08817

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DGUOK

Uniprot

Q16854

Synonyms:

deoxyguanosine kinase; DGUOK; dGK ;

Immunogen:

Human DGUOK.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A