

DHCR24 Antibody



PACO08818

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DHCR24

Uniprot

Q15392

Synonyms:

24-dehydrocholesterol reductase; DHCR24; KIAA0018; Nbla03646; SELADIN1; seladin-1 ;

Immunogen:

Human DHCR24.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A