

DHRS1 Antibody



PACO08823

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DHRS1

Uniprot

Q96LJ7

Synonyms:

dehydrogenase/reductase (SDR family) member 1; DHRS1; DKFZp586I0523; FLJ14250; FLJ25430; MGC20204; SDR19C1 ;

Immunogen:

Human DHRS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A