

# DHX32 Antibody



PACO08832

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DHX32

**Uniprot**

Q7L7V1

**Synonyms:**

DEAH (Asp-Glu-Ala-His) box polypeptide 32; DHX32; DDX32; DHLP1; FLJ10694; FLJ10889 ;

**Immunogen:**

Human DHX32.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A