

DHX9 Antibody



PACO08836

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DHX9

Uniprot

Q08211

Synonyms:

DEAH (Asp-Glu-Ala-His) box polypeptide 9; DHX9; DDX9; LKP; NDHII; RHA ;

Immunogen:

Human DHX9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A