

DIRAS2 Antibody



PACO08842

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

DIRAS2

Uniprot

Q96HU8

Synonyms:

DIRAS family, GTP-binding RAS-like 2; DIRAS2; DKFZp761C07121; Di-Ras2 ;

Immunogen:

Human DIRAS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A