

DLD Antibody



PACO08850

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DLD

Uniprot

P09622

Synonyms:

dihydrolipoamide dehydrogenase; DLD; DLDH; E3; GCSL; LAD; PHE3 ;

Immunogen:

Human DLD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A