

DLEU1 Antibody



PACO08852

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DLEU1

Uniprot

O43261

Synonyms:

deleted in lymphocytic leukemia 1 (non-protein coding); DLEU1; BCMS; DLB1; DLEU2; LEU1; LEU2; MGC22430; NCRNA00021; XTP6 ;

Immunogen:

Human DLEU1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A