

DLX4 Antibody



PACO08860

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DLX4

Uniprot

Q92988

Synonyms:

distal-less homeobox 4; DLX4; BP1; DLX7; DLX8; DLX9 ;

Immunogen:

Human DLX4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A