

DNAI1 Antibody



PACO08868

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNAI1

Uniprot

Q9UI46

Synonyms:

dynein, axonemal, intermediate chain 1; DNAI1; CILD1; ICS; ICS1; MGC26204; PCD ;

Immunogen:

Human DNAI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A