

DNAI2 Antibody



PACO08869

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DNAI2

Uniprot

Q9GZS0

Synonyms:

dynein, axonemal, intermediate chain 2; DNAI2; CILD9 ;

Immunogen:

Human DNAI2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A