

DNAJA1 Antibody



PACO08870

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DNAJA1

Uniprot

P31689

Synonyms:

DnaJ (Hsp40) homolog, subfamily A, member 1; DNAJA1; DJ-2; DjA1; HDJ2; HSDJ; HSJ2; HSPF4; hDJ-2 ;

Immunogen:

Human DNAJA1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A