

DNAJB11 Antibody



PACO08875

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DNAJB11

Uniprot

Q9UBS4

Synonyms:

DnaJ (Hsp40) homolog, subfamily B, member 11; DNAJB11; ABBP-2; ABBP2; DJ9; EDJ; ERdj3; ERj3; HEDJ; PRO1080; UNQ537; hDj9 ;

Immunogen:

Human DNAJB11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A