

DNAJB8 Antibody



PACO08883

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNAJB8

Uniprot

Q8NHS0

Synonyms:

DnaJ (Hsp40) homolog, subfamily B, member 8; DNAJB8; DJ6; MGC33884 ;

Immunogen:

Human DNAJB8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A