

DNAJC1 Antibody



PACO08884

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DNAJC1

Uniprot

Q96KC8

Synonyms:

DnaJ (Hsp40) homolog, subfamily C, member 1; DNAJC1; DNAJL1; ERdj1; HTJ1; MGC131954; MTJ1 ;

Immunogen:

Human DNAJC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A