

DNAJC12 Antibody



PACO08886

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DNAJC12

Uniprot

Q9UKB3

Synonyms:

DnaJ (Hsp40) homolog, subfamily C, member 12; DNAJC12; JDP1 ;

Immunogen:

Human DNAJC12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A